

Rendr Portal Roles 2.4 - External and Internal

Action \ Role	EXTERNAL					INTERNAL		
	Facility Admin	Review Doctor	Lead Tech	Field Tech	Office Pers	Support	Production	Super Admin
Dashboard	●	●	●	●	●	●		●
View Patient	●	S	●	S	●	●		●
Create Patient	●		●	●	●	●		●
Edit Patient	●		●	●	●	●		●
Add Documents	●	S	●	●	●	●		●
View Studies List	●	S	●	S	●	●		●
View (Open) Study	●	S	●	S		●		●
Edit Study	●	S	●	S		●		●
Share Study	●	S	●	S		●		●
View Shared Users	●	S	●	●		●		●
View Studies Audit Log	●		●			●		●
Mark Study for Deletion	●		●			●		●
Import Study	●		●	●		●		●
View User	●		●			●		●
Invite User	●					●		●
Create User	●					●		●
View User Audit Logs	●		●			●		●
Edit User	●					●		●
Delete User	●					●		●
Create Amplifier	CONTACT LIFELINES SUPPORT TO ADD AN ACQUISITION DEVICE TO YOUR FACILITY					●		●
View Amplifier	●					●	●	●
Create Device	CONTACT LIFELINES SUPPORT TO ADD AN ACQUISITION DEVICE TO YOUR FACILITY					●	●	●
View Device						●	●	●
Create Cameras	CONTACT LIFELINES SUPPORT TO ADD AN ACQUISITION DEVICE TO YOUR FACILITY					●	●	●
View Cameras	●					●	●	●
View Multiple Facilities*	●*	●*	●*	●*	●*	●	●	●
Edit Facility	●							●
Montages	●	●	●	●			●	
Help	●	●	●	●	●	●	●	●

* - USER MUST BE A USER IN THE FACILITIES
 S - Shared or Created - the user can do this if the study is shared with him/her or if they created the study.
 ● - The user has full permission for the function.