

F is for Flexible EEG Solution

Flexibility for In-home or Physician Office EEG

Outpatient EEG, Streaming Review with Rendr Platform Cloud Software

The Trackit F includes user-friendly features designed with input from neurodiagnostic clinicians. Remotely monitor patients with streaming data on multiple web browsers with Rendr Platform cloud software. Accessible on multiple devices including PC, Mac and tablet for study access anywhere, anytime.



Outpatient Ambulatory Video EEG

The Trackit F EEG System is a complete solution for outpatient EEG and HD video monitoring. While it can be used for ambulatory studies in both clinical and home settings, it was designed for a workflow in which the EEG tech deploys the test in a patient's home. The durable plastic hardshell case with customized moldings is rugged and sanitizable, created for in-home use.

Features Include:

- 22" x 14" x 9" high impact protection travel case with handles and wheels
- Waterproof seal (when closed)
- Case holds all components, including amplifier, photic, camera and card reader
- Durable closed-cell case for easy cleaning and protection from scent absorption
- Ultra HD PTZ (pan, tilt, zoom) video camera with 120° clamp mount, with remote control via Rendr software to ensure patients are kept on camera




Rendr Platform Cloud Software

Streaming data review, from anywhere on multiple devices

The Trackit F EEG System streams data and HD video to the Rendr Platform via the cloud, with a robust internet connection. Physicians can review live studies from anywhere, on any device, manipulating montage and filters independently from the technologist. This allows immediate feedback and rapid evaluation of potential seizure activity.

Data Protection

Amazon Web Services (AWS), the host of the Rendr Platform, includes HIPAA-ready architecture. AWS includes load sensing and auto-scaling to meet increased bandwidth demand so that as your patient load grows your ability to record and review their studies will too.



EEG Anywhere Subscription

Our flexible subscription plan features the Trackit F bundled with our Rendr Platform software for one fixed monthly fee. With unlimited seat licensing and rapid replacement if your system needs repair, our EEG Anywhere Subscription provides predictable monthly billing, with no large up front capital expenditure, to rapidly scale your business. As your partner, we'll work with you to scale and expand the number of patients you serve.

